

# A sterile, single use filter



## Sterifilt®

**Packaging:** Individual sterile blister pack.

**Composition:** Polypropylene membrane filter.

**Properties:** Reduces injection related harms, such as: abscess, phlebitis, “dirty hits”, embolism and oedema. Removes 90% of insoluble particles bigger than 5 microns, 99% of particles 10 microns and larger. Retains virtually no drug residue.

**Indications:** Prevention of tissue damage and infections related to the unintentional injection of insoluble particles. Promotes the single-use of filters. Helps reduce risk of blood-borne virus transmission.

**Caution:** Sterifilt degrades and may become blocked after use. Re-use should not be attempted. Sterifilt does not filter out viruses or bacteria.

**Clinical history:** Developed in France since 1997 in conjunction with needle exchange programs. Safety and efficacy validated by laboratory assays. Tested with brown heroin (base), white heroin (hydrochloride), cocaine (hydrochloride), crack (cocaine base), amphetamines and crushed tablets.

**Directions for use:** Designed for use with 1 ml Insulin syringes (Becton Dickinson, BBraun, Nevershare Syringe, Terumo, Artsana, Chirana) or any type of syringe (1 ml, 2 ml, 5 ml, 10 ml).

If you are using a syringe that doesn't have a needle attached, the Sterifilt fits where the needle goes. Fix the needle once you have drawn up the drug and removed the used sterifilt. Hold vertically. Draw up the liquid slowly. Safely dispose of filter and syringe immediately after use, using approved sharps-safe or similar puncture and shatter-proof container.

## Sterifilt® is compatible with all types of syringes

Sterifilt is a CE registered medical device.

Sterifilt



# Order Form

Complete this form and post to:

APOTHICOM Distribution/52 avenue Edison. 75013 Paris. FRANCE  
Tél.: 00 33 (0)1 53 61 18 41. Fax: 00 33 (0)1 53 61 04 49. E-mail: [infos@apothicom.org](mailto:infos@apothicom.org)

Please fill in using capital letters

NAME: ..... FIRST NAME: .....

JOB TITLE: ..... ORGANISATION: .....

DELIVERY ADDRESS: .....

.....

INVOICE ADDRESS (if different): .....

.....

PHONE NUMBER: ..... FAX NUMBER: .....

E-MAIL ADDRESS: .....@ .....

To make ordering simpler, Sterifilt® is available and priced, in multiples of 1,000. Sterifilt® is either displayed in an individual sterile blister or inserted inside the Stericup® blister set instead of the cotton filter. This latter version, called "Sterimix®", has been designed to offer a more complete and effective sterile tool.

We aim to deliver promptly; if you require an urgent delivery please call. Orders can be placed by fax to: +33 1 53 61 04 49 or transmitted by mail through a PDF file at [infos@apothicom.org](mailto:infos@apothicom.org)



Sterifilt® (individual blister) Shipping costs are not included				Quantity	Total
From 1 000 to 5 000	Per thousand	Ex Vat 120,00	Tax incl. 143,52		
6 000 and more	Per thousand	Ex Vat 115,00	Tax incl. 137,54		
Sterifilt Leaflets	10 per 1 000	Free			Free
Sterifilt Poster	1 per 1 000	Free			Free
Sterifilt information cards	20 per 1 000	Free			Free
Shipping costs are not included					
Sterimix® (Stericup® + Sterifilt® + dry swab + handle)				Quantity	Total
From 1 000 to 5 000	Per thousand	Ex Vat 195,00	Tax incl. 233,22		
6 000 and more	Per thousand	Ex Vat 190,00	Tax incl. 227,24		
Sterifilt Leaflets	10 per 1 000	Free			Free
Sterifilt Poster	1 per 1 000	Free			Free
Sterifilt information cards	20 per 1 000	Free			Free
Shipping costs are not included					

All items are priced ex shipping. Transport costs can be quoted according to quantities, please let us know if you would need a prior quotation. For the European Union countries only: if you are VAT registered please specify your intracommunautary number on the above form, so that we shall not add French VAT rate (19,6%) to the invoice amount.

Payment Method: by transfer bank or by cheque to Apothicom.

Date:

Signature:

Update notice: May 2009. Please cancel all previous and invalid notices

For further information • e-mail: [infos@apothicom.org](mailto:infos@apothicom.org) or call + 33 (0)1 53 61 18 41

**Bank**

Name: Apothicom Distribution  
IBAN Number: FR76 3000 4003 1700 0100 6193 565  
BIC: BNPAFRPPBY